

Stay in the
HOME
you love!



An independent
family business
since 1867



Stannah
The Stairlift People

Stannah – a name you can *trust*

Stannah is a family-owned business with over 140 years of experience in the lifting industry. We are proud to make some of the world's best-selling stairlifts, with over 500,000 sold globally. We welcome the opportunity to share that experience with you, either in person or through our network of the best stairlift dealers in North America. A stairlift is an important purchase, so it makes sense to choose Stannah – a name you can trust.

Better by design

A good product depends partly on using quality materials, but far more on the skill of designers and engineers. Stannah has a worldwide reputation for thoughtful design; our stairlifts are styled to fit into your home, are easy to use, reliable, and have appropriate safety features included as standard.

Examples of our thoughtful design include the way our stairlifts can be folded without bending down (especially important at the top of the stairs), the seat belt placement on the end of an armrest (much easier to reach than a seat belt by the hips) and the shape of the track which is narrow, and enclosed so that anything dropped on it falls harmlessly aside.

All of our stairlifts use a reliable battery system with a smart trickle-charger that plugs into a standard power outlet. Not only does this use very little electricity, it is also quiet and makes for a smooth ride that will allow you to use your stairlift in the event of a power outage.

Safety

We believe that safety is the most important factor when making a stairlift. That's why we include more standard safety features than any other stairlift manufacturer and insist on high standards of work when we are training installers. In addition, every stairlift we design is tested for the equivalent of ten years of use and the motor is tested for 100,000 journeys (which is about 20 years of normal residential use) to ensure that it is reliable enough to meet your needs.

However, engineering is only part of the story when it comes to safety. Just as important are the needs of the user and the staircase in question. That's why it is so important to have an expert come and review your situation before buying a stairlift.



Room to Maneuver

The 600 model will fit your straight stairs perfectly, leaving plenty of room when the lift isn't in use.

Shown above: Siena chair, Tan upholstery.



Enjoy the Great Outdoors!

The Stannah 320 outdoor model is designed to withstand harsh weather, salty spray and small animals.



Stairs That Turn?

At Stannah, we'll custom design our Model 260 to fit your stairs. The rail can also be set to finish in the most convenient position.

Shown above: Starla chair, Stone upholstery.

Product Features



Seat Belt

We offer four seat belt options, including a retractable, secure belt that fastens and unfastens easily with one hand.



Safety Sensors

Every Stannah stairlift is equipped with sensors that detect obstructions and bring your unit to a safe and immediate stop.



Swivel Options

Whether you choose the power option or the standard manual swivel, your lift will turn to let you off in the safest position.



Hand Control

Easy to use, the slimline hand control allows you to use any part of your hand to operate your stairlift.



Armrest Control

The conveniently placed On-Off switch is clearly marked and the power footrest button (where fitted) makes folding your chair safe and easy.



Remote Control

Each stairlift comes standard with two remote controls, allowing you to call your chair to the top or bottom landing.



Easy Folding

All of our indoor stairlifts come with footrest folding assistance, so you'll find that your lift is both safe and unobtrusive.



Upholstery Choices

We offer several upholstery colors and styles, including wood trim options for Starla chairs.



Powered Folding Rail

This option for straight stairs folds up the rail automatically, leaving plenty of walking space on the lower landing.



“How nice it was to find a company more interested in what’s best for me and not just best for them”

Stannah. One simple call is all it takes.

Contact your local Stannah supplier, who will arrange to visit your home for a friendly, free consultation. They will listen to your needs, answer your questions and provide a full, no-obligation quotation. If you decide to go ahead, they will measure your staircase using the latest technology to style your stairlift.

Stannah promises to provide:

- Your rail, custom-made to fit your stairs
- Your stairlift, tailored to your requirements
- Quick and clean installation
- Expert support, close at hand
- Comprehensive, peace-of-mind guarantee

Your Stannah supplier.

Call for a free no-obligation survey of your stairs.



6300 Falls Road - Unit 2 Baltimore, MD 21209

Tel: (410) 825-1440

9210 Baltimore Ave. College Park, MD 20740

Tel: (301) 585-0700

Email: Info@bedcomobility.com

Our company policy is one of continuing product improvement. We therefore reserve the right to alter specifications and improve design at any time. Stannah is a registered trademark. Colors may vary due to normal printing processes.

©Stannah Stairlifts Inc 2014. Not to be reproduced wholly or in part without the written permission of Stannah Stairlifts Inc.

Stannah

Stannah Stairlifts Inc.

Head office & showroom:

20 Liberty Way, Suite A
Franklin MA 02038

Toll-free: 1-800-877-8247

Tel: 1-508-520-7878

Fax: 1-508-520-7883

E-mail: Stannah@Stannah.com

New York showroom:

45 Knollwood Rd., Suite 103
Elmsford, NY 10523

Toll-free: 1-888-333-8708

Tel: 1-914-345-8484

Fax: 1-914-345-8488

E-mail: Stannah@Stannah.com

www.Stannah-Stairlifts.com

MA HIC #160211. CT Elevator Limited Contractor License ELV0475333-R5. Not to be reproduced wholly or in part without the written permission of Stannah Stairlifts Inc.