




Ricon Wheelchair Lifts™



Clearway™

Power.....electro-hydraulic	Rated load capacity800 lbs*
Motor rating @ 12 volts DC 65.0 amp avg/cycle, 1250 psi	Manual backup (up) hand pump
Hydraulic cylinders 2 ea, 1.5", power-up/gravity-down	Manual backup (down).....pressure release valve
	Lift weight 340 - 370 lbs



Specifications (inches)

	A	B	C	D	E	F	G	H
MODEL	Stationary frame width	Height (folded)	Installation depth (folded)	Usable platform width	Usable platform length	Floor-to-ground travel	Traveling frame width	Clear entry width
S1132	41	40	19.5	26	38	27	33.5	27
S1133	45	40	21.5	30	38	27	37.5	31
S1230	45	50	21.5	30	48	31	37.5	31
S1231	45	44	21.5	30	42	31	37.5	31

* Lift is rated at 800 lbs on Dodge Promaster and Ford Transit when installed with vehicle specific install kits.