



## TECHNICAL SPECIFICATIONS

**Lift Capacity:** 350 lb (159 kg)

**Platform Width:** 28.5 in (72.4 cm)

**Max. Mobility Device Width:** 28 in (71 cm outside of wheel to outside of wheel)

**Hitch class:** II or III\*

**Vehicle Type:** CUVs, hatchbacks, minivans, pickups, sedans, SUVs, full-size vans

*\*Receiver hitch required*

Learn more at [www.bruno.com](http://www.bruno.com)

Keep all of your vehicle's interior space for seating and cargo when you transport your scooter or powerchair with Bruno's Out-Sider exterior lift. It's durable, safe and reliable. Easily gain access to cargo space by adding the "Swing-Away" option. Available for nearly every type of vehicle.

- Fully-powered, single-button operation
- Automatic folding platform when not in use
- No scooter or powerchair modifications
- Low profile design for rear view mirror visibility
- Hold-Tite foot or retractable self-tensioning belts for mobility device securement
- Integrated safety latch to prevent unintentional lift lowering
- Manual backup for peace of mind
- Potential to reinstall in another applicable vehicle\*
- Eligible for auto manufacturer potential rebates\*
- 3-year limited warranty

*\* Contact Bruno dealer for details.*

## OUT-SIDER A/B FEATURES



Hold-Tite foot securement for scooters (ASL-250A).



Retractable belt securement for powerchairs/scooters (ASL-250B).



Self-leveling platform with traction coating.



One-button operation with built-in safety latch to prevent unintentional platform lowering.

## OUT-SIDER A/B OPTIONS

- Vinyl cover (not shown)



Low-profile third brake light.



Swing-Away: Easy access to hatch/cargo space.



Bruno Independent Living Aids, Inc.  
1780 Executive Drive, Oconomowoc, WI 53066  
[www.bruno.com](http://www.bruno.com)

© Bruno Independent Living Aids, Inc. 2019  
2019-04-10-Out-Sider-A-B.indd • MKT-L-3